



FEATURES

- PMMA Electrolytic Cell (volume 256ml)
- Integated Rectifier 10A -12V (resolution: 0,01A - 0,1V), Ripple 1%, V-A stabilization
- Electric Heater to heat up the liquid
- Adjustable Mechanical Agitation
- Adjustable Insufflation System
- Integrated Timer with acoustic alarm

APPLICATIONS

Analysis of the quality of the plating bath:

- Adhesion of the coating to the metal
- Deposit Internal Tensions
- Metal contamination
- Additives Balancing and quantity
- Current density Optimization

ADVANTAGES

- All the instruments in one Structure
- PPS flameproof construction
- Removable cell-cart with drip proof system



SPECIFICATION

VERSION LEGEND			
	Insufflation	Agitator	Heater
001	✓	✗	✗
002	✓	✓	✓
003	✓	✓	✗




TESTER IDENTIFICATION					
Model	Rectifiers N°	Voltage (Volt)	Frequency (Hz)	Material	Version (see the legend)
GTM	R1=1	220	50 60	PS=PPS	001
	R2=2				002
					003
GTM	R1	220	50	PS=PPS	002


DIMENSIONS



GTC CELL



Cell
GTC SZ001



Insufflation Cell
GTC IN001



CELL IDENTIFICATION

Description	Type	Version
GTC	IN=with insufflation	001
	SZ=without insufflation	
GTC	IN	001

Product Range: PUMP: Magnetic Drive, Horizontal, Mechanical Seal, Vertical, Diaphragm, Drum, - FILTERING SYSTEMS – PRECIOUS METAL RECOVERY SYSTEMS – DEPURATION – FILTERS: String Wound Cartridges, Melt Blown, Filtering Disks, Bag Filters, Plated Filter Cartridges, Absolute – ACCESSORIES: Gratings, Venturi Eductor – CUSTOMIZED COMPONENTS.

LAFONTE.EU may modify data any time, without notice, in order to get product enhancement